



# MILK DOCTOR BETA

## Screening test kits for Beta Lactam Tests In Bovine Milk

- ▶ Unnecessary extraction the samples
- ▶ Specifically match to Japanese MRL level
- ▶ Result in 5 minutes
- ▶ All materials in one kit
- ▶ Simple operation and user friendly



MADE  
IN  
JAPAN

50test

# SIMPLE!

# EASY!

# FASTNESS!

**USER FRIENDLY**

Unnecessary to extract samples

**FASTNESS**

Result in 5 minutes

**NO INSTRUMENTS**

Unnecessary incubator to heat samples

**HIGH SENSITIVITY**

Penicillin 3ppb  
Ceftiofur 5ppb  
Cloxacillin 5ppb

**EYES INTERPRETATION**

Easy to get result

## Typical antibiotics Detection Limits

- Penicillin	3ppb
- Cefazolin	3ppb
- Nafcillin	5ppb
- Cefroxime	5ppb
- Dicloxacillin	20ppb
- Cefutiofur	5ppb
- Cefalonium	3ppb
- Cephapirin	20ppb
- Cloxacillin	4ppb
- Cefquinom	20ppb



## Example of result


Negative

Positive



## Operating Procedure

Collect milk samples



Add milk samples



Dip test paper and wait at room temperature



AR Brown Co., Ltd Specialty Chemical dept.

Agri Food safety FD-B team

Dai-ni Marutaka Bldg., 7-13-8, Ginza, Chuo-ku, Tokyo, Japan

TEL : 81-3-3545-5724 FAX : 81-3-3543-8865

URL : <http://www.arbrown.com/> MAIL: [a-r@arbrown.com](mailto:a-r@arbrown.com)